

Cleaning, Care & Maintenance



IMPORTANT CLEANING INFORMATION

- Use only 3% Hydrogen Peroxide Solution for routine cleaning and disinfection.
- If hydrogen peroxide is unavailable, use warm water and a mild dish detergent only (no additives, fragrances, lotions, antibacterial agents, bleach, or harsh chemicals).
- Do NOT use heavy detergents, alcohol, bleach, solvents, or abrasive cleaning tools.
- When cleaning the system, do not wet the outside of the Collection Bag Hydro-Block™ Air Filter.
- Always allow components to air dry completely before storage.

DAILY SYSTEM CLEANING **URINCARE CUP™ SYSTEM**

Step 1 – Flush the System

Leave all components connected:

- Collection Bag
- Control Unit
- URINCARE CUP™
 1. Fill the **CUP** halfway with 3% hydrogen peroxide solution.
 2. The **Control Unit** may activate automatically. If it does not:
 - Press and hold the blue "M" button until all liquid has been evacuated.
 3. Repeat this process twice.
 4. Repeat the same procedure using clean water to rinse the system.
 5. Empty the **Collection Bag** completely.
 6. Ensure the drain valve is fully closed when finished.



Component Cleaning

Step 2 – Clean Individual Components

Disconnect the System

1. Disconnect the **CUP** from the **Collection Bag** hose by squeezing the connector release.
2. Press the black release button to remove the **Control Unit** from the **Collection Bag** clip.

Clean the Control Unit

1. Fill a small container with 3% hydrogen peroxide solution.
2. Submerge the first ½ inch of the connector end of the **Control Unit**.
3. Press the blue "M" button and run the pump for approximately 4 seconds.
4. Repeat using clean water.

Clean the CUP

1. Place the connector end of the **CUP** into 3% hydrogen peroxide solution for approximately 4 seconds.
2. Remove and allow to air dry.
3. Repeat using clean water if desired.

Cleaning Urine Contamination

Step 3 – Control Unit & Battery Cleaning

If urine contacts the **Control Unit** or **Rechargeable Battery Pack**:

- Press the battery release and remove the **battery** pack.
- Using a paper towel or soft cloth dampened with 3% hydrogen peroxide, wipe all contaminated surfaces.
- Inspect **battery** contacts and connector pins.
- If corrosion is present:
 - Use a soft toothbrush lightly dampened with 3% hydrogen peroxide.
 - Gently brush the **battery** contacts and mating contacts on the **Control Unit**.
- Wipe away any residue with a clean paper towel.
- Allow all components to dry completely before reconnecting.

Cleaning, Care & Maintenance



Collection Bag Care

- Empty the Collection Bag regularly.
- Inspect tubing and connectors for buildup or blockage.
- Verify the Hydro-Block™ Air Filter remains dry and unobstructed.
- Replace Collection Bags and accessories according to your prescribed replacement schedule.

Routine Inspection

Before each use, inspect:

- Tubing for kinks, tears, or blockage.
- Connectors for secure attachment.
- Battery contacts for cleanliness and corrosion.

Collection Bag and drain valve for leaks.

Discontinue use of damaged components and contact Omni Medical Systems for replacement parts.

Questions or Support?

For replacement supplies, technical support, or product assistance, please contact Omni Medical Systems Customer Support at 1-802-497-2253 (Option 1) or customercare@urincare.com.



Connector End of Male Cup that attaches to the Control Unit

Where Male Cup connects to the Control Unit on the Collection Bag



**Male Harness
Machine Washable**



**Hydro-Block™
Air Filter**

**Collection Bag
Drain Valve**