

HARNESS SIZE	WAIST SIZE
XS / SMALL	24" - 30"
MED / LARGE	32" - 38"
X-LARGE	40" - 42"
2X / 3X	44" - 50"



***The URINCARE™ System is not indicated for patients with a penile length of less than 1.5 inches or for patients with a retracted penis.*

PROPER POSITIONING AND HARNESS ADJUSTMENT

For best results, the **URINCARE CUP™** should be held firmly against the body in a near-vertical position. Proper positioning is essential for effective urine collection and device performance.

The harness features adjustable **Velcro®** straps that can be tightened or loosened to achieve a secure, comfortable fit. A properly adjusted harness helps maintain the correct position and seal of the cup.

“Tuck the Cup”: Always ensure the penis is pointed downward and fully tucked into the CUP before use and the CUP is pointed down.

Important: The CUP should remain firmly seated against the body throughout use to ensure reliable operation, as the sensor are located at the bottom of the CUP.

